

Product Description

LoFric Elle, a single use hydrophilic urinary catheter, designed for female users.

Sterilisation and shelf life

Sterilised by irradiation. Use before expiry date on the package.

Packaging

30 catheters per customer box. 2 customer boxes per shipping box.

Material composition

Catheter Material:
Polyolefin-based elastomer (POBE).

Surface coating:
Polyvinyl pyrrolidone (PVP).

Wetting fluid:
Water-based wetting fluid with glycerol.

Primary package Material:
Polypropylene.


Product range

Reference	Type	Length	Size
418 10 25	LoFric® Elle™ Nelaton (female)	10 cm*	CH 10
418 12 25	LoFric® Elle™ Nelaton (female)	10 cm*	CH 12
418 14 25	LoFric® Elle™ Nelaton (female)	10 cm*	CH 14

*Effective length, length of the catheter that can be inserted into the body.



User-friendly customer box.



- Designed to be ergonomic, user-friendly, with an angulated handle – for good grip and non-touch technique. Connects to a urine collection bag or extension tube.
- Ready to use – Ready when you need it, to enable fast handling and catheterisation.
- Slim and Discreet design.
- Lower container is recyclable as plastic (polypropylene).
- Not manufactured with PVC, Phthalates and latex.
- Smooth eyelets.^{1,2}
- LoFric catheters are coated with the unique Urotonic™ Surface Technology. This ensures the catheter remains smooth and safe, both during insertion and withdrawal.^{3,4}

1. Lundgren et al. 2000
2. Newman et al. 2018

3. Waller et al 1997
4. Dunger et al. 1988

At Wellspect we develop innovative continence care solutions that change people's lives. We are committed to inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 30 years with our product brands LoFric® and Navina™. We create reliable and user-friendly products for bladder and bowel management with as little environmental impact as possible. We passionately strive to become climate neutral and work closely together with users and healthcare professionals who constantly inspire us to improve our products and services in a sustainable way, now and for the future.

Wellspect. A real difference.

For more information about our products and our initiative Advancing Continence Care Together (ACCT), please visit [Wellspect.com](https://www.wellspect.com).

Join the conversation on Facebook and Instagram.

[wellspect.co.uk](https://www.wellspect.co.uk)



Rx Only