



Loop opening.

Adhesive tab on back for reliable hanging.

Instant activation - sterile salt solution included, squeeze and it's ready to use.

Protective sleeve for optimal control allows for non-touch technique.

Not manufactured with PVC, phthalates or latex.



LoFric catheters are coated with the unique Urotonic™ Surface Technology. This ensures the catheter remains smooth and safe, both during insertion and withdrawal.<sup>1,2</sup>



Foldable into handy, discreet pocket size.

### Product description

LoFric® Origo™, single use, hydrophilic urinary catheter with salt solution for activation.

### Sterilisation and shelf life

Sterilised by irradiation. Use before expiry date on package.

### Package

30 catheters per customer box.  
4 customer boxes per shipping box.

5 catheters per sample pack.

### Material composition

**Hydrophilic surface coating:** Polyvinylpyrrolidone (PVP).

**Catheter:** Polyolefin-based elastomer (POBE).

**Sacht:** Laminate foil of polyethylene terephthalate (PET), aluminum and polyethylene (PE). Sterile water with sodium chloride (NaCl).

**Sleeve:** Polyethylene (PE)

### Item packaging:

Top foil: Laminate foil of polyethylene terephthalate (PET) and polyethylene (PE).

Bottom foil: Polyolefin.

### Product range

Reference	Type	Length	Size
446 12 25	LoFric® Origo™ Sleeve Nelaton (male)	40 cm	CH 12
446 14 25	LoFric® Origo™ Sleeve Nelaton (male)	40 cm	CH 14
446 16 25	LoFric® Origo™ Sleeve Nelaton (male)	40 cm	CH 16

All products available in sample pack, 5 catheters.  
To order replace the last two digits in the Reference number with 05.



User-friendly packaging.



Sample pack, 5 catheters.



**LoFric® Origo™**

## References

1. Waller et al., Spinal Cord, 1997
2. Dunger et al., Acta Endocrinologica, 1988
3. Lundgren et al., Spinal Cord, 2000
4. Newman et al., Springer, 2018

At Wellspect we develop innovative continence care solutions that improve quality of life for people with bladder and bowel problems. We inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 40 years with our product brands LoFric® and Navina™. We always aim to minimise the environmental impact of our products and passionately strive to become climate neutral. We work together with users and healthcare professionals to improve clinical outcome in a sustainable way, now and for the future.

**Wellspect. A Real Difference.**

For more information about our products and services, please visit [Wellspect.co.uk](https://www.wellspect.co.uk).

Join the conversation on LinkedIn, Twitter, Facebook and Instagram.

[wellspect.co.uk](https://www.wellspect.co.uk)

Wellspect Limited, Brunel Way, Stonehouse, Gloucestershire GL10 3GB.  
Phone: 0800 652 3350. [www.wellspect.co.uk](https://www.wellspect.co.uk)

LoFric



Manufacturer: Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.  
Phone: +46 31 376 40 00. [www.wellspect.com](https://www.wellspect.com)

CE 2797

Rx Only